

Recombinant Klebsiella pneumoniae Kanamycin kinase type II,NEO

Catalog No: #AP76268

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP76268-1 10ug #AP76268-2 50ug #AP76268-3 500ug #AP76268-4 1mg

Description

Product Name	Recombinant Klebsiella pneumoniae Kanamycin kinase type II,NEO
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Klebsiella pneumoniae
Accession No.	P00552
Uniprot	P00552
GeneID	39540766;
Calculated MW	29kDa
Target Sequence	MIEQDGLHAGSPAAWVERLFGYDWAQQTIGCSDAAVFRLSAQGRPVLVKTDLSGALNELQDEAARLSWLAT TGVPCAAVLDVVTEAGRDWLLLGEVPGQDLLSSHLAPAEKVSIMADAMRRLHTLDPATCPFQAKHRIERA RTRMEAGLVDQDDLDEEHQGLAPAEFLARLKARMPDGEDLVVTHGDACLPNIMVENGRFSGFIDCGRLGVAD RYQDIALATRDIAEELGGEWADRFLVLYGIAAPDSQRIAFYRLLDEFF
Formulation	Supplied as a 0.2 µm filtered solution of PBS, pH 7.4, 20% Glycerol.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only