

Recombinant Human Interleukin-22,IL-22 (C-hlgG2 Fc)

Catalog No: #AP76269

Orders: order@signalwayantibody.com

Package Size: #AP76269-1 10ug #AP76269-2 50ug #AP76269-3 500ug #AP76269-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Interleukin-22,IL-22 (C-hlgG2 Fc)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9GZX6
Uniprot	Q9GZX6
GeneID	50616;
Calculated MW	43.4kDa
Target Sequence	<p>APISSHCRLDKSNFQQPYITNRTFMLAKEASLADNNTDVRLIGEKLFGVSMSERCYLMKQVLNFTLEEVLPQ</p> <p>SDRFQPYMQEVVFLARLSNRLSTCHIEGDDLHIQRNVQKLKDTVKKLGESGEIKAIGELDLLFMSLRNACIGG</p> <p>GGSGGGGGSGGGGSERKCCVECPCPAPPVAGPSVFLFPPKPKDTLMISRTPEVTCVVVDVSHEDPEVQFNW</p> <p>YVDGMEVHNAKTKPREEQFNSTFRVSVLTVVHQDWLNGKEYKCKVSNKGLPAPIEKTISKTKGQPREPQVY</p> <p>TLPPSREEMTKNQVSLTCLVKGFYPSDIAVEWESNGQPENNYKTTTPMLDSDGSFFLYSKLTVDKSRWQQG</p> <p>NVFSCVMHEALHNHYTQKSLSLSPGK</p>
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only