

Recombinant Mouse OX40 Ligand, TNFSF4, OX40L (N-8His)

Catalog No: #AP76273

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP76273-1 10ug #AP76273-2 50ug #AP76273-3 500ug #AP76273-4 1mg

Description

Product Name	Recombinant Mouse OX40 Ligand, TNFSF4, OX40L (N-8His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	P43488
Uniprot	P43488
GeneID	22164;
Calculated MW	17.7kDa
Target Sequence	HHHHHHHHSSSPAKDPPIQRLRGAVTRCEDGQLFISSYKNEYQTMEVQNNSVVIKCDGLYIIYLKGSFFQEVKI DLHFREDHNPISIPMLNDGRRIVFTVVASLAFKDKVYLVNAPDTLCEHLQINDGELIVQLTPGYCAPEGSYHS TVNQVPL
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only