

Recombinant E. coli Tryptophan Synthase (N-6His)

Catalog No: #AP76274



Package Size: #AP76274-1 10ug #AP76274-2 50ug #AP76274-3 500ug #AP76274-4 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant E. coli Tryptophan Synthase (N-6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	E.coli,E.coli
Accession No.	P0A877&P0A879
Uniprot	P0A877
GeneID	946204;
Calculated MW	72.5kDa
Target Sequence	MERYESLFAQLKERKEGAFVFPVTLGDPGIEQSLKIIDLIEAGADALELGIPFSDPLADGPTIQNATLRAFAAGV TPAQCFEMLALIRQKHPTIPIGLLLMYANLVFNKGIDEFYQCEKVGVDVSVLVADVPVEESAPFRQAALRHNVAPI FICPPNADDDLLRQIASYGRGYTYLLSRAGVTGAENRAALPLNHLVAKLKEYNAAPPLQGFGISAPDQVKAID AGAAGAISGSAIVKIIIEQHINEPEKMLAALKVFVQPMKAATRS&MHHHHHHTLLNPNYFGEFGGMYVPQILMPA LRQLEEFVSAQKDPEFQAQFNDLLKNYAGRPTALTKCQNITAGTNTTLYLKREDLLHGGAHKTNQVLGQALL AKRMGKTEIIAETGAGQHGVASALASALLGLKCRIYMGAKDVERQSPNVFRMLMGAEVIPVHSGSATLKDAC NEALRDWSGSYETAHYMLGTAAGPHPYPTIVREFQRMIGEETKAQILEREGRLPDAVIACVGGGSNAIGMFAD FINETNVGLIGVEPGGHGIETGEHGAPLKHGRVGIYFGMKAPMMQTEDGQIEESYSISAGLDFPSVGPQHAYL NSTGRADYVSITDDEALEAFKTLCLHEGIIPALESSHALAHALKMMRENPDKEQLLVNLSGRGDKDIFTVHDIL KARGEI
Formulation	Supplied as a 0.2 µm filtered solution of PBS, pH7.4.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only