

Recombinant Human Interleukin-36γ,IL-36γ,IL-1F9

Catalog No: #AP76275



Package Size: #AP76275-1 10ug #AP76275-2 50ug #AP76275-3 500ug #AP76275-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Interleukin-36γ,IL-36γ,IL-1F9
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9NZH8
Uniprot	Q9NZH8
GeneID	56300;
Calculated MW	17kDa
Target Sequence	SMCKPITGTINDLNQQVWTLQGQNLVAVPRSDSVTPVTVAVITCKYPEALEQGRGDPIYLGIQNPMECLYCEK VGEQPTLQLKEQKIMDLYGQPEPVKPFIFYRAKTGRTSTLESVAFPDWFIASSKRDQPIILTSELGKSYNTAFEL NIND
Formulation	Lyophilized from a 0.2 μm filtered solution of 20mM Tris,100mM Nacl,0.1mM EDTA,pH8.0.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only