

# Recombinant Human C-C Motif Chemokine 2,CCL2,MCP-1



Catalog No: #AP76276

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP76276-1 10ug #AP76276-2 50ug #AP76276-3 500ug #AP76276-4 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human C-C Motif Chemokine 2,CCL2,MCP-1
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P13500
Uniprot	P13500
GeneID	6347;
Calculated MW	8.7kDa
Target Sequence	QPDAINAPVTCCYNFTNRKISVQRLASYRRITSSKCPKEAVIFKTIVAKEICADPKQKWVQDSMDHLDKQTQTP KT
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only