

Recombinant Human Cadherin-1,E-Cadherin,CDH1 (C-6His)



Catalog No: #AP76283

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP76283-1 10ug #AP76283-2 50ug #AP76283-3 500ug #AP76283-4 1mg

Description

Product Name	Recombinant Human Cadherin-1,E-Cadherin,CDH1 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P12830
Uniprot	P12830
GeneID	999;
Calculated MW	61.2kDa
Target Sequence	DWVIPPISCPENEKGFPPKNLVQIKSNKDKEGKVFYSITGQGADTPPVGVFIERETGWLKVTPLDRERATYT LFSHAVSSNGNAVEDPMEILITVTDQNDNKPEFTQEVFKGSVMEGALPGTSVMEVTATDADDDVNTYNAAIAY TILSQDPELDPKMNFTINRNTGVISVTTGLDRESFPTYTLVVQAADLQGEGLSTTATAVITVTDNDNPIFNPT TYKGQVPENEANVVITTLKVTADAPNTPAWEAVYTILNDDGGQFVVTNPNVNDGILKTAKGLDFEAKQQYIL HVAVTNVVPFEVSLTTSTATVTDVLDVNEAPIFVPEKRVESDFGVGQEITSYTAQEPDTFMEQKITYRIW RDTANWLEINPDTGAISTRAELDREDFEHVKNSTYTALIIATDNGSPVATGTGTTTTLSDVNDNAIPEPRTIFFC ERNPKPQVINIIDADLPPNTSPFTAELTHGASANWTIQYNDPTQESIILKPKMALEVGDYKINLKLMDNQNKDQV TTLEVSVCDCEGAAGVCRKAQPVEAGLQIHVVVVVV
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only