

Recombinant Human IL-23 Recetor,IL-23R (C-6His)

Catalog No: #AP76309



Package Size: #AP76309-1 10ug #AP76309-2 50ug #AP76309-3 500ug #AP76309-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human IL-23 Recetor,IL-23R (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q5VWK5
Uniprot	Q5VWK5
GeneID	149233;
Calculated MW	38.7kDa
Target Sequence	GITNINCSGHIWVEPATIFKMGMNISYQAAIKNCQPRKLHFYKNGIKERFQITRINKTTARLWYKNFLEPHASM YCTAECPKHFQETLICGKDISSGYPPDIPDEVTCVIYEYSGNMTCTWNAGKLTIDTKYVVHVKSLETEEEQQY LTSSYINISTDSLQGGKKYLVWVQAANALGMEESKQLQIHLDDIVIPSAAVISRAETINATVPKTIYWDSQTTIEK VSCEMRYKATTNQTWNVKEFDTNFTYVQQSEFYLEPNIKYVFQVRCQETGKRYWQPWSSPFFHKTPETVPQ VTSKAFQHDTWNSGLTVASISTGHILTSDNRGDHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only