

# Recombinant Human Complement Component C3a,C3a



Catalog No: #AP76310

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP76310-1 10ug #AP76310-2 50ug #AP76310-3 500ug #AP76310-4 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human Complement Component C3a,C3a
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P01024
Uniprot	P01024
GeneID	718;
Calculated MW	9.1kDa
Target Sequence	SVQLTEKRMDKVGKYPKELRKCCEDGMRENPMRFSCQRRTRFISLGEACKKVFLDCCNYITELRRQHARASH LGLAR
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only