

Recombinant Macaca mulatta (Rhesus macaque) α-2-HS-Glycoprotein,AHSG,Fetuin A (C-10His)

Catalog No: #AP76313

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP76313-1 10ug #AP76313-2 50ug #AP76313-3 500ug #AP76313-4 1mg

Description

Product Name	Recombinant Macaca mulatta (Rhesus macaque) α-2-HS-Glycoprotein,AHSG,Fetuin A (C-10His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Macaca mulatta,Rhesus Macaque
Accession No.	F6WRS6
Uniprot	F6WRS6
Calculated MW	38.9kDa
Target Sequence	APHGPGLIYRQPNCDDEPETEEAALVAVDYINQNLPGWYKHTLNQIDEVKVWPQQPFGEIIEIDTLETTCHV LDPTPVAGCRVRLKEHAVEGDCDFQLLKVNGKFSVEYAKCDSSPDSAEDVRKVCRCDFLLAPLNDTRVVHA AEAAVTAFNAQNNGSNFQLEEISRAQLVPLPPSTYVEFTISGTDCAKEATEAANCNLLAKKQYGFCKATLNEK LGGEVAVTCTVFQTQPATSQPQPEGANEAVPTPVVDADAPASYPVGASGPLPTGSPYAHVLAAPPVPIH SSHYDLRHLFMGVVSLRSPSGEALHPRKTRSVVQPSVGAAGPVVPPCPGRVRHFKVHHHHHHHHHH
Formulation	Lyophilized from a 0.2 μm filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only