

Recombinant Human Complement Factor D,Adipsin (C-10His)

Catalog No: #AP76314

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP76314-1 10ug #AP76314-2 50ug #AP76314-3 500ug #AP76314-4 1mg

Description

Product Name	Recombinant Human Complement Factor D,Adipsin (C-10His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P00746
Uniprot	P00746
GeneID	1675;
Calculated MW	25.8kDa
Target Sequence	ILGGREAEAHARPYMASVQLNGAHLCCGGVLVAEQWVLSAAHCLEDAADGKVQVLLGAHSLSQPEPSKRLYD VLRAVPHPDSPDPTIDHDLQLLSEKATLGPAVRPLPWQRVDRDVPAGTLCDVAGWGIVNHAGRRPDSLQH VLLPVLDRATCNRRTHHDGAITERLMCAESNRRDCKGDSGGPLVCGGVLEGVVTSGSRVCGNRKKPGIYTR VASAAWIDSVLAHHHHHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM Tris, 150mM NaCl, pH7.5.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only