

Recombinant Human IL-1 Receptor Type 2,IL-1R-2 (C-Fc)



Catalog No: #AP76316

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP76316-1 10ug #AP76316-2 50ug #AP76316-3 500ug #AP76316-4 1mg

Description

Product Name	Recombinant Human IL-1 Receptor Type 2,IL-1R-2 (C-Fc)
Host Species	Human cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P27930
Uniprot	P27930
GeneID	7850;
Calculated MW	64.5kDa
Target Sequence	FTLQPAHTGAARSCRFRGRHYKREFRLEGEFVALRCPQVPYWLWASVSPRINLTWHKNSARTVPGEEET RMWAQD GALWLLPALQEDSGTYVCTTRNASYCDKMSIELRVFENTDAFLPFISYPQILTLSTSGVLVCPDLSEF TRDKTDVKIQWYKDSLLLDKDNEKFLSVRGTHLLVHDVALEDAGYYRCVLTFAHEGQQYNITRSIELRIK EETIPVIISPLKTISASLGSRLTIPCKVFLGTGTPLTMLWWTANDTHIESAYPGGRVTEGPRQEYSENNENYIEV PLIFDPVTRELDHMDFKCVVHNTLSFQTLRRTVKEIEGRDMDPKSCDKTHTCPPCPAPELLGGPSVFLFPPKPK DTLMISRTPEVTCVVVDVSHEDPEVKFNWYVDGVEVHNAKTKPREEQYNSTYRVVSVLTVLHQDWLNGKEYK CKVSNKALPAPIEKTISKAKGQPREPQVYTLPPSRDELTKNQVSLTCLVKGFYPSDIAVEWESNGQPENNYKT TPPVLDSDGSFFLYSKLTVDKSRWQQGNV FSCSVMHEALHNHYTQKSLSLSPGK
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only