

Recombinant Mouse Cytotoxic T-Lymphocyte Protein 4,CTLA-4,CD152 (C-6His)



Catalog No: #AP76319

Orders: order@signalwayantibody.com

Package Size: #AP76319-1 10ug #AP76319-2 50ug #AP76319-3 500ug #AP76319-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Mouse Cytotoxic T-Lymphocyte Protein 4,CTLA-4,CD152 (C-6His)
Host Species	Human cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	P09793
Uniprot	P09793
GeneID	12477;
Calculated MW	14.6kDa
Target Sequence	AIQVTQPSVVLASSHGVASFPCEYSPSHNTDEV RVTVLRQTNDQMTEVCATTFTEKNTVGFLDYPFCSGTFNE SRVNLTIQGLRAVD TGLYLCKVELMYPPYFVGMGNGTQIYVIDPEPCPDS DHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only