

Recombinant Mouse IL-1 Receptor Type 1,IL-1R-1 (C-Fc)



Catalog No: #AP76320

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP76320-1 10ug #AP76320-2 50ug #AP76320-3 500ug #AP76320-4 1mg

Description

Product Name	Recombinant Mouse IL-1 Receptor Type 1,IL-1R-1 (C-Fc)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	P13504
Uniprot	P13504
GeneID	16177;
Calculated MW	64kDa
Target Sequence	LEIDVCTEYPNQIVLFLSVNEIDIRKCLTPNKMHGDTIIWYKNSKTPISADRDSRIHQNEHLWFVPAKVEDS GYYYCIVRNSTYCLKTKVTVTLENDPGLCYSTQATFPQRLHIAGDGLVCPYVSYFKDENNELPEVQWYKNC KPLLLDNVSFFGVKDKLLVRNVAEEHRGDYICRMSYTFRGKQYPVTRVIQFITIDENKRDRPVILSPRNETIEAD PGSMIQLICNVTGQFSDLVYWKWNGSEIEWNDPFLAEDYQFVEHPSTKRKYTLITLNLISEVKSQFYRYPFICV VKNTNIFESAHVQLIYPVPDFKIEGRDMDPKSCDKTHTCPPCPAPELLGGPSVFLFPPKPKDTLMISRTPEVTC VVVDVSHEDPEVKFNWYVDGVEVHNAKTKPREEQYNSTYRVVSVLTVLHQDWLNGKEYKCKVSNKALPAPIE KTISKAKGQPREPQVYTLPPSRDELTKNQVSLTCLVKGFYPSDIAVEWESNGQPENNYKTPPVLDSDGSFFL YSKLTVDKSRWQQGNVFCFSVMHEALHNHYTQKSLSLSPGK
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only