

Recombinant Mouse IL-1 Receptor Type 2,IL-1R-2 (C-Fc)

Catalog No: #AP76321

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP76321-1 10ug #AP76321-2 50ug #AP76321-3 500ug #AP76321-4 1mg

Description

Product Name	Recombinant Mouse IL-1 Receptor Type 2,IL-1R-2 (C-Fc)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	P27931
Uniprot	P27931
GeneID	16178;
Calculated MW	64.9kDa
Target Sequence	FTTPTVVHTGKVSESPITSEKPTVHGDNCQFRGREFKSELRLEGEVVLRCPLAPHSDISSSSHSFLTWSKLD SSQLIPRDEPRMWWKGNILWILPAVQQDSGTICTFRNASHCEQMSVELKVKNTEASLPHVSYLQISALSTTG LLVCPDLKEFISSNADGKIQWYKGAILLDKGNKEFLSAGDPTRLLISNTSMDDAGYYRCVMTFTYNGQEYNITR NIELRVKGTTEPIPIVISPLETIPASLGSRLIVPCKVFLGTGTSSNTIVWWLANSTFISAAAYPRGRVTEGLHHQYS ENDENYVEVSLIFDPVTREDLHTDFKCVASNPRSSQLHTTVKEIEGRDMDPKSCDKTHTCPPCPAPELLGGP SVFLFPPKPKDTLMISRTPEVTCVVVDVSHEDPEVKFNWYVDGVEVHNAKTKPREEQYNSTYRVVSVLTVLH QDWLNGKEYKCKVSNKALPAPIEKTISKAKGQPREPQVYTLPPSRDELTKNQVSLTCLVKGFYPSDIAVEWES NGQPENNYKTTPVLDSDGSFFLYSKLTVDKSRWQQGNVFCFSVMHEALHNHYTQKSLSLSPGK
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only