

## Recombinant Rhesus B7-2,CD86 (C-6His)

Catalog No: #AP76323



Package Size: #AP76323-1 10ug #AP76323-2 50ug #AP76323-3 500ug #AP76323-4 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Rhesus B7-2,CD86 (C-6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Cynomolgus monkey, Cynomolgus monkey
Accession No.	H9ZF18
Uniprot	H9ZF18
Calculated MW	25.7kDa
Target Sequence	YFNETADLPCQFANSQNRSLSELVFWQEQENLVLNEVYLGKEKFDSVHSKYMGRTSFDPESWTLRLHNLQI KDKGLYQCIIHHRPTGMIRIHQMNSELSVLNFSQPEIVPISNITENMYINLTCSSIHGYPEPEKMSVLLRTKNS TIEYDGVMMQKSQDNVTELYDVSISLSVSFPDVTSNMTIFCVLETDKTQLLSSPFSIELEDPQPPPDHIPPHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only