

Recombinant Rabbit B7-2,CD86 (C-6His)

Catalog No: #AP76324



Package Size: #AP76324-1 10ug #AP76324-2 50ug #AP76324-3 500ug #AP76324-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Rabbit B7-2,CD86 (C-6His)
Host Species	Human cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Rabbit
Accession No.	P42071
Uniprot	P42071
GeneID	100008999;
Calculated MW	26.3kDa
Target Sequence	AASLRIQAYFNKTADLPCQFTNSQSRSLSELVVFWDQERLVLVLYELFLGREKPDNVDPKYIGRTSFDQESWNL QLHNVQIKDKGVYQCFVHHRGAKGLVPIYQMNSELSVLNFTQPEITLISNITRNSAINLTCSSVQGYPEPKMF FVLKTENATTEYDGVIEKSQDNVTGLYNISISGSITFSDDIRNATIYCVLQTESTETYSQHFPIVPADPVPVEKPR HHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only