

Recombinant Cynomolgus B7-1,CD80 (C-6His)

Catalog No: #AP76336



Package Size: #AP76336-1 10ug #AP76336-2 50ug #AP76336-3 500ug #AP76336-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Cynomolgus B7-1,CD80 (C-6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Cynomolgus monkey, Cynomolgus monkey
Accession No.	G7MKE2
Uniprot	G7MKE2
Calculated MW	24.7kDa
Target Sequence	VIHVTKEVKEVATLSCGHNVSEELAQTRIYWQKEKKMVLTMMSGDMNIWPEYKNRTIFDITNNLSIVILALRPS DEGTYECVVLKYEKDAFKREHLAEVMLSVKADFPTPSITDFEIPPSNIRRIICSTSGGFPEPHLSWLENGEELNA ISTTVSQDPETELYTVSSKLDNFNMTTNHSFMCLIKYGHRLRVNQTFNWNTPKQEHFPDNHHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only