

Rad51 (Phospho-Thr309) Polyclonal Antibody

Catalog No: #12383



Package Size: #12383-1 50ul #12383-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rad51 (Phospho-Thr309) Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms
Specificity	Phospho-Rad51 (T309) Polyclonal Antibody detects endogenous levels of Rad51 protein only when phosphorylated at T309.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human Rad51 around the phosphorylation site of T309.
Target Name	Rad51
Modification	Phospho
Other Names	RAD51; RAD51A; RECA; DNA repair protein RAD51 homolog 1; HsRAD51; hRAD51; RAD51 homolog A
Accession No.	Swiss-Prot: Q06609NCBI Gene ID: 5888
Uniprot	Q06609
GeneID	5888;
Target Species	human
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1:500-2000

Immunohistochemistry: 1/100 - 1/300

ELISA: 1/5000

Images



Western blot analysis of HELA-UV 293T-UV using p-Rad51 (T309) antibody.

Note: This product is for in vitro research use only