

Recombinant Mouse B7-1,CD80 (C-6His)

Catalog No: #AP76337



Package Size: #AP76337-1 10ug #AP76337-2 50ug #AP76337-3 500ug #AP76337-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Mouse B7-1,CD80 (C-6His)
Host Species	Human cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	Q00609
Uniprot	Q00609
GeneID	12519;
Calculated MW	24.6kDa
Target Sequence	VDEQLSKSVKDKVLLPCRYNSPHEDESEDRIYWQKHDKVLSVIAGKLVWPEYKNRTLYDNNTYSLIILGLVL SDRGTYSCVVQKKERGTYEVKHLALVKLSIKADFSTPNITESGNPSADTKRITCFASGGFPKPRFSWLENGRE LPGINTTISQDPESELYTISSQLDFNTRNHTIKCLIKYGDAHVSEDFTWEKPPEDPPDSKNHHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only