

Recombinant Rat T-lymphocyte Activation Antigen CD80,B7-1 (C-6His)

Catalog No: #AP76338

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP76338-1 10ug #AP76338-2 50ug #AP76338-3 500ug #AP76338-4 1mg

Description

Product Name	Recombinant Rat T-lymphocyte Activation Antigen CD80,B7-1 (C-6His)
Host Species	Human cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Rat
Accession No.	O35187
Uniprot	O35187
Calculated MW	25.7kDa
Target Sequence	VGLFQISSGIVGQVSKSVREKALLSCDYKFCSEEQSIHRIYWQKHDKMVLVSVISGVPEVWPEYKNRTVYDIANN YFSLLGLILSDRGTYTCVVQRYEGESYVVKHLTTVELSVRADFPTPNITESGNPSADIKRITCFASGGFPKPRL SWLENGRELNGINNTISQDPESELYTISSQLDFNTTYDHFIDCFIEYGDAHVSQNFTWVKPPEDPPDEKQH ^{HHH} HH
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only