

## Recombinant Rabbit B7-1,CD80 (C-6His)

Catalog No: #AP76339



Package Size: #AP76339-1 10ug #AP76339-2 50ug #AP76339-3 500ug #AP76339-4 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Rabbit B7-1,CD80 (C-6His)
Host Species	Human cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Rabbit
Accession No.	P42070
Uniprot	P42070
GeneID	100009377;
Calculated MW	24.5kDa
Target Sequence	GISQVTKSVKEMAALSCDYNISIDELARMRIYWQKDQQMVLISIISGQVEVWPEYKNRTFPDIINNLSLMILALRLS DKGTYTCVVQKNENGSRREHLTSVTLRIRADFPVPSITDIGHPDPNVKRIRCSASGGFPEPRLAWMEDGEEL NAVNTTVDQDLDELYSVSELDNFVNTNHSIVCLIKYGELSVSQIFPWSKPKQEPPIDQHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only