

Recombinant Human Importin Subunit Beta-1,KPNB1(N-6His)



Catalog No: #AP76340

Orders: order@signalwayantibody.com

Package Size: #AP76340-1 10ug #AP76340-2 50ug #AP76340-3 500ug #AP76340-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Importin Subunit Beta-1,KPNB1(N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q14974
Uniprot	Q14974
GeneID	3837;
Calculated MW	98.6kDa
Target Sequence	MNHKVVHHHHHMHMELITILEKTVSPDRLELEAAQKFLERAAVENLPTFLVELSRVLANPGNSQVARVAAGLQIKN SLTSKDPDIKAQYQQRWLAIDANARREVKNYVLQTLGTETYRPSSASQCVAGIACAEIPVNWPELIPQLVANV TNPNSTEHMKESTLEAIGYICQDIDPEQLQDKSNEILTAIQGMRKEEPSNNVLAATNALLNSLEFTKANFDKE SERHFIMQVVCEATQCPDTRVRVAALQNLVKIMSLYYQYMETYMGPALFAITIEAMKSDIDEVALQGIEFWSNV CDEEMDLAIEASEAAEQGRPPEHTSKFYAKGALQYLVPILTQTLTKQDENDDDDWNPCAAGVCLMMLLATC CEDDIVPHVLPFIKEHIKNDWRYRDAAVMAFGCILEGPEPSQLKPLVIQAMPTLIELMKDPSVVVRDTAAWTV GRICELLPEAAINDVYLAPLLQCLIEGLSAEPRVASNVCFWAFSSLAEEAAYEADVADDQEEPATYCLSSSFELIV QKLETTDRPDGHQNNLRSSAYESLMEIVKNSAKDCYPAVQKTTLVIMERLQQVLQMESHQSTSDRIQFNDL QSLLCATLQNVLRKQHQDALQISDVVMASLLRMFQSTAGSGGVQEDALMAVSTLVEVLGGEFLKYMEAFKP FLGIGLKNYAEYQVCLAAGVGLVGDLCRALQSNIPFCDEVMLLLENLGNENVHRSVKPQILSVFGDIALAIGGE FKKYLEVVLNLTQQASQAQVDKSDYDMVDYLNELRESCLEAYTGIVQGLKGDQENVHPDVMLVQPRVEFILSF IDHIAGDEDHTDGVVACAAGLIGDLCTAFGKDVLLKLEARPMIHELLTEGRRSKTNKAKTLATWATKELRKLKN QA
Formulation	Supplied as a 0.2 µm filtered solution of 20mM Tris,1mM DTT,30% glycerol,0.1M Nacl,pH8.0.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only