

Recombinant Mouse T-cell Surface Protein Tactile,CD96(C-6His)



Catalog No: #AP76342

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP76342-1 10ug #AP76342-2 50ug #AP76342-3 500ug #AP76342-4 1mg

Description

Product Name	Recombinant Mouse T-cell Surface Protein Tactile,CD96(C-6His)
Host Species	Human cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	Q3U0X8
Uniprot	Q3U0X8
GeneID	84544;
Calculated MW	57.6kDa
Target Sequence	VWEELFNVGDDVYALPGSDINLTCQTKEKNFLVQMOWSKVTDKNDMIALYHPQYGLYCGQEHACESQVAAT ETEKGVTNWTLYLRNISSALGGKYECIFTLYPEGIKTTVYNLIVEPYTQDEHNYTIEIETNRTLEIPCFQNTSSEIP PRFTFSWLVEKDGVEEVLFTHHHHVNNSTSFKGRIRLGGDYRLHLSVPVQIQDDGRTF SCHLTVNPLKAWKMS TTVKVFAKPEILMTVENSTMDVLGERVFTCLLKNVFPKANITWFIDGRFLQGNEEGIYITNEEKNCSSGFWELK SVLTRMHSGPSQSNMATAWCMALSPGPRNKMWNTSSQPITVSFDSVIAPTKHLPTVTGSTLGTQPFSDAGVS PTGYLATPSVTIVDENGLTPDATPQTSNSSMTTKDGNYLEASSGTDAKNSSRAAASSKSGSWPFPFTSPEW HSLPGTSTGPEPDPSPVSWIPSEVHTSAPLDASLAPHDTIISTTTEFPNVLTTANGTTKIDHGPITSIIVNQPSDG MHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only