

TIRAP (Phospho-Tyr86) Polyclonal Antibody

Catalog No: #12384



Package Size: #12384-1 50ul #12384-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TIRAP (Phospho-Tyr86) Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms
Specificity	Phospho-TIRAP (Y86) Polyclonal Antibody detects endogenous levels of TIRAP protein only when phosphorylated at Y86.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human TIRAP around the phosphorylation site of Y86.
Target Name	TIRAP
Modification	Phospho
Other Names	TIRAP; MAL; Toll/interleukin-1 receptor domain-containing adapter protein; TIR domain-containing adapter protein; Adaptor protein Wyatt; MyD88 adapter-like protein
Accession No.	Swiss-Prot: P58753NCBI Gene ID: 114609
Uniprot	P58753
GeneID	114609;
Target Species	human
Calculated MW	23kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1:500-2000

Immunohistochemistry: 1/100 - 1/300

ELISA: 1/5000

Images

Western blot analysis of K562 using p-TIRAP (Y86) antibody.



Note: This product is for in vitro research use only