

Recombinant Cynomolgus Tumor Necrosis Factor Ligand 2B,OX40L(N-8His)



Catalog No: #AP76347

Orders: order@signalwayantibody.com

Package Size: #AP76347-1 10ug #AP76347-2 50ug #AP76347-3 500ug #AP76347-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Cynomolgus Tumor Necrosis Factor Ligand 2B,OX40L(N-8His)
Host Species	Human cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Cynomolgus monkey
Accession No.	F7FL80
Uniprot	F7FL80
Calculated MW	16.5kDa
Target Sequence	HHHHHHHHQVSHQYPRIQSIKVQFTEYKKEEGFILTSQKEDEIMKVQNNNSVIINCDGFYLISLKGYFSQEVNISL HYQKDEEPLFQLKKVRSVNSLMVASLTYKDKVYLVNTTDDNTSLDDFHVNGGELILIHQNPGEFCVL
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only