

TBC1D4 (Phospho-Thr642) Polyclonal Antibody

Catalog No: #12385

Package Size: #12385-1 50ul #12385-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TBC1D4 (Phospho-Thr642) Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Hu Ms
Specificity	Phospho-TBC1D4 (T642) Polyclonal Antibody detects endogenous levels of TBC1D4 protein only when phosphorylated at T642.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human TBC1D4 around the phosphorylation site of T642.
Target Name	TBC1D4
Modification	Phospho
Other Names	TBC1D4; AS160; KIAA0603; TBC1 domain family member 4; Akt substrate of 160 kDa; AS160
Accession No.	Swiss-Prot: O60343NCBI Gene ID: 9882
Uniprot	O60343
GeneID	9882;
Target Species	human
SDS-PAGE MW	150kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1:500-2000

Immunohistochemistry: 1/100 - 1/300

Immunofluorescence: 1/200 - 1/1000

ELISA: 1/5000

Images



Western blot analysis of L929 using p-TBC1D4 (T642) antibody.

Note: This product is for in vitro research use only