

## Recombinant Mouse TRANCE,RANK L,TNFSF11

Catalog No: #AP76365



Package Size: #AP76365-1 10ug #AP76365-2 50ug #AP76365-3 500ug #AP76365-4 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Mouse TRANCE,RANK L,TNFSF11
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	O35235
Uniprot	O35235
GeneID	21943;
Calculated MW	20.1kDa
Target Sequence	QRFSGAPAMMEGSWLDVAQRGKPEAQPFHAHLTINAASIPSGSHKVTLSWYHDRGWAKISNMTLSNGKLRV NQDGFYLYLANICFRHHETSGSVPTDYLQLMVYVVKTSIKIPSSHNLKGGSTKNWGSNSEFHFYSINVGFF KLRAGEEISIQVSNPSSLDPDQDATYFGAFKVQDID
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only