

Recombinant Human Fibroblast Growth Factor 4,FGF-4

Catalog No: #AP76367

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP76367-1 10ug #AP76367-2 50ug #AP76367-3 500ug #AP76367-4 1mg

Description

Product Name	Recombinant Human Fibroblast Growth Factor 4,FGF-4
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P08620
Uniprot	P08620
GeneID	2249;
Calculated MW	16.9kDa
Target Sequence	MSLARLPVAAQPKEAAVQSGAGDYLLGIKRLRRLYCNVGIGFHLQALPDGRIGGAHADTRDSLLELSPVERGV VSIFGVASRFFVAMSSKGLYGSPFFTDECTFKEILLPNYNAYESYKYPGMFIALSKNGKTKKGNRVSPMTMKV THFLPRL
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only