

Recombinant Mouse C-C Motif Chemokine 3,CCL3,MIP-1a(N-6His)

Catalog No: #AP76370

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP76370-1 10ug #AP76370-2 50ug #AP76370-3 500ug #AP76370-4 1mg

Description

Product Name	Recombinant Mouse C-C Motif Chemokine 3,CCL3,MIP-1a(N-6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	P10855
Uniprot	P10855
GeneID	20302;
Calculated MW	10.8kDa
Target Sequence	MGSSHHHHHHSSGLVPRGSHMDDDDKAPYGADTPTACCFYSYRKIPRQFIVDYFETSSLCSQPGVIFLTARN RQICADSKETWVQEYITDLELNA
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM Tris, 150mM NaCl, 5% Trehalose, 1mM EDTA, pH8.0.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only