

Recombinant Human Artemin

Catalog No: #AP76373



Package Size: #AP76373-1 10ug #AP76373-2 50ug #AP76373-3 500ug #AP76373-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Artemin
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q5T4W7
Uniprot	Q5T4W7
GeneID	9048;
Calculated MW	12.1kDa
Target Sequence	GSAGGPGSRARAAGARGCRLRSQ LVPVRALGLGHRSDLVRFRCGSGCRRARSPHDLASLLGAGALRP PPGSRPVSPCCRPTRYEA VSFMDVNSTWRTVDRLSATA CGCLG*
Formulation	Lyophilized from a 0.2 µm filtered solution of 4mM HCl.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only