

## Recombinant Mouse Galectin-9 (N-GST)

Catalog No: #AP76378



Package Size: #AP76378-1 10ug #AP76378-2 50ug #AP76378-3 500ug #AP76378-4 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Mouse Galectin-9 (N-GST)
Host Species	E. coli
Purification	Greater than 90% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	O08573-2
Uniprot	O08573
GeneID	16859;
Calculated MW	62.8kDa
Target Sequence	MSPILGYWKIKGLVQPTRLLEYLEEKYEEHLYERDEGDKWRNKKFELGLEFPNLPYYIDGDVKLTQSMAIIRYI ADKHNMLGGCPKERAEISMLEGAVLDIRYGVSR IAYSKDFETLKVDFLSKLP EMLKMFEDRLCHKTYLNGDHV THPDFMLYDALDVVLYMDPMCLDAFPKLVCFKKRIEAIQIDKYLKSSKYIAWPLQGWQATFGGGDHPPKSDL VPRGSMALFSAQSPYINPIIPFTGPIQGGLEGLQVTLQGTTKSFAQRVVFQNSFNNGNDIAFHFNPRFEEGG YVVCNTKQNGQWGPEERKMMPFQKMPFELCFLVQRSEFKVMVNKKFFVQYQHRVPYHLVDTIAVSGCL KLSFITFQTQNFRAHQAPMAQTIIHMHVHSTPGQMFSTPGIPPVVYPTPAYTIPFYTPINGLYPSKSIMISGNV LPDATRFHINLRCGGDIAFHNLNPRFNENAVVRNTQINNSWGQEERSLLGRMPFSRGQSFSVWIIICEGHCFKVA VNGQHMCEYYHRLKNLQDINTLEVAGDIQLTHVQT
Formulation	Supplied as a 0.2 µm filtered solution of 1M PBS, 15% glycerol,pH8.0.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only