

# Recombinant Mouse ICOS Ligand,B7-H2,CD275 (C-6His)

Catalog No: #AP76383

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP76383-1 10ug #AP76383-2 50ug #AP76383-3 500ug #AP76383-4 1mg

## Description

Product Name	Recombinant Mouse ICOS Ligand,B7-H2,CD275 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	Q9JHJ8
Uniprot	Q9JHJ8
GeneID	50723;
Calculated MW	27.1kDa
Target Sequence	ETEVGAMVGSNNVLSCIDPHRRHFNLGSLYVYWQIENPEVSVTYLPLYPKSPGINVDSSYKNRGHLSLDSMKQ GNFSLYLKNVTPQDTQEFTCRVFMNTATELVKILEEVRLRVAANFSTPVISTDSSNPGQERTYTCMSKNGY PEPNLYWINTTDNSLIDTALQNNTVYLNKLGLYDVISTLRPLPWTSRGLVCCVENVALHQNITSISQAESFTGNN TKNPQETHNNELKHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only