

Recombinant Human Fc γ RIIIA,FCGR3A,CD16a (C-6His,Val176Phe)

Catalog No: #AP76386

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP76386-1 10ug #AP76386-2 50ug #AP76386-3 500ug #AP76386-4 1mg

Description

Product Name	Recombinant Human Fc γ RIIIA,FCGR3A,CD16a (C-6His,Val176Phe)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P08637
Uniprot	P08637
GeneID	2214;
Calculated MW	22.7kDa
Target Sequence	GMRTEDLPKAVVFLEPQWYRVLEKDSVTLKCGAYSPEDNSTQWFHNESLISSQASSYFIDAATVDDSGEYR CQTNLSTLSDPVQLEVHIGWLLLQAPRWVFKEEDPIHLRCHSWKNTALHKVTYLQNGKGRKYFHHNSDFYIPK ATLKDSGSYFCRGLFGSKNVSSETVNITITQGLAVSTISSFFPPGYQHSHHHH
Formulation	Lyophilized from a 0.2 μ m filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only