

Recombinant Mouse Fc γ RIIIA,FCGR3,CD16 (C-6His)

Catalog No: #AP76400

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP76400-1 10ug #AP76400-2 50ug #AP76400-3 500ug #AP76400-4 1mg

Description

Product Name	Recombinant Mouse Fc γ RIIIA,FCGR3,CD16 (C-6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	P08508
Uniprot	P08508
Calculated MW	22kDa
Target Sequence	ALPKAVVKLDPPWIQVLKEDMVTLMCEGTHNPGNSSTQWFHNGRSIRSQQASYTFKATVNDSGEYRCQME QTRLSDPVDLGVISDWLLLQTPQRVFLEGETITLRCHSWRNKLLNRISFFHNEKSVRYHHYKSNFSIPKANHSH SGDYCYCKGSLGSTQHQSKPVTITVQDPATTSSISLVWYHTHHHHHH
Formulation	Lyophilized from a 0.2 μ m filtered solution of PBS, pH7.5.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 4 months.

Note: This product is for in vitro research use only