

Recombinant Human Interleukin-4 Receptor Subunit Alpha,IL-4 R α (C-6His)



Catalog No: #AP76404

Orders: order@signalwayantibody.com

Package Size: #AP76404-1 10ug #AP76404-2 50ug #AP76404-3 500ug #AP76404-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Interleukin-4 Receptor Subunit Alpha,IL-4 R α (C-6His)
Host Species	Human cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P24394
Uniprot	P24394
GeneID	3566;
Calculated MW	24.4kDa
Target Sequence	MKVLQEPTCVSDYMSISTCEWKMNPTNCSTELRLLYQLVFLLEAHTCIPENGGAGCVCHLLMDDVVSAD NYTLDLWAGQQLLWKSFKPSEHVKPRAPGNLTVHTNVSDTLTLLTWSNPYPDPNLYNHLTYAVNIWSENDP ADFRIYNVTYLEPSLRIAASTLKSGISYRVRVRAWAQCYNNTWSEWSPSTKWHNSYREPFEQHSHHHHH
Formulation	Lyophilized from a 0.2 μ m filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only