

Recombinant Human B7 Homolog 6,B7-H6,NCR3LG1(C-6His)

Catalog No: #AP76410

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP76410-1 10ug #AP76410-2 50ug #AP76410-3 500ug #AP76410-4 1mg

Description

Product Name	Recombinant Human B7 Homolog 6,B7-H6,NCR3LG1(C-6His)
Host Species	Human cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q68D85
Uniprot	Q68D85
GeneID	374383;
Calculated MW	27.5kDa
Target Sequence	DLKVEMMAGGTQITPLNDNVTIFCNIFYSQLNITSMGITWFWKSLTFDKEVKVFEFFGDHQEAFRPGAIVSPW RLKSGDASLRLPGIQLLEEAGEYRCEVVVTPKLAQGTQVQLEVVASPARLLLDQVGMKENEDKYMCESSGFYP EAINITWEKQTQKFPHPHIESEDVITGPTIKNMDGTFNVTSCCLKLNSSQEDPGTVYQCVRHASLHTPLRSNFTL TAARHSLSETEKTDNFSHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only