

Recombinant Cynomolgus TIM-3,HAVCR2 (C-Fc)

Catalog No: #AP76421



Package Size: #AP76421-1 10ug #AP76421-2 50ug #AP76421-3 500ug #AP76421-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Cynomolgus TIM-3,HAVCR2 (C-Fc)
Host Species	Mammalian cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Cynomolgus monkey, Cynomolgus monkey
Accession No.	G7P6Q7
Uniprot	G7P6Q7
Calculated MW	46.3kDa
Target Sequence	SEVEYIAEVGQNAAYLPCSYTPAPPGNLVPVCWVGKACPVFDCSNVLRDNRDNDRTSGRYWLKGFHKG DVSLTIENVTLADSGVYCCRIQIPGIMNDEKHNVKLVVIKPAKVTPAPTLQRDLTSAFPRMLTTGEHGAETQTP GSLPDVNLTVSNFFCELQIFLTNLRD SGATIRIEGRMDPKSCDKTHTCPPCPAPELLGGPSVFLFPPKPKDT LMISRTPEVTCVWVDVSHEDPEVKFNWYVDGVEVHNAKTKPREEQYNSTYRVVSVLTVLHQDWLNGKEYKC KVSNAKALPAIEKTISKAKGQPREPQVYTLPPSRDELTKNQVSLTCLVKGFYPSDIAVEWESNGQPENNYKTP PVLDSGDGSFFLYSKLTVDKSRWQQGNVFSCSVMHEALHNHYTQKSLSLSPGK
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only