

Recombinant Human Mucin-1, MUC-1 (C-Fc)

Catalog No: #AP76425



Package Size: #AP76425-1 10ug #AP76425-2 50ug #AP76425-3 500ug #AP76425-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Mucin-1, MUC-1 (C-Fc)
Host Species	Human Cells
Purification	Greater than 90% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P15941-11
Uniprot	P15941
GeneID	4582;
Calculated MW	42.3kDa
Target Sequence	APKPATVVTGSGHASSTPGGEKETSATQRSSVPSSTEKNAFNSSLEDPSTDYYQELQRDISEMFLQIYKQGGF LGLSNIKFRPGSVVVQLTLAFREGTINVHDVETQFNQYKTEAASRYNLTISDVSVDVPPFFSAQSGAGVPGDI EGRMDPKSCDKTHTCPPCPAPELLGGPSVFLFPPKPKDTLMISRTPEVTCVVDVSHEDPEVKFNWYVDGVE VHNAKTKPREEQYNSTYRVVSVLTVLHQDWLNGKEYKCKVSNKALPAPIEKTISKAKGQPREPQVYTLPPSRD ELTKNQVSLTCLVKGFYPSDIAVEWESNGQPENNYKTTTPVLDSDGSFFLYSKLTVDKSRWQQGNVFCFSVM HEALHNHYTQKSLSLSPGK
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only