

# Recombinant Mouse Thrombopoietin, THPO, TPO (C-6His)

Catalog No: #AP76426

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP76426-1 10ug #AP76426-2 50ug #AP76426-3 500ug #AP76426-4 1mg

## Description

Product Name	Recombinant Mouse Thrombopoietin, THPO, TPO (C-6His)
Host Species	Human Cells
Purification	Greater than 90% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	P40226
Uniprot	P40226
GeneID	21832;
Calculated MW	36.4kDa
Target Sequence	SPVAPACDPRLLNKLLRDSHLLHSRLSQCPDVPDPLSIPVLLPAVDFSLGEWKTQTEQSKAQDILGAVSLLLEGV MAARGQLEPSCLSLLGQLSGQVRLLLGALQGLLGTQLPLQGRRTTAHKDPNALFLSLQQLLRGKVRFLLLVEG PTLCVRRRLPTTAVPSSTSQLLTLNKFPNRTSGLLETNFSVTARTAGPGLLSRLQGFRVKITPGQLNQTSRSPV QISGYLNRTHGPNVNGTHGLFAGTSLQTLASDISPGAFNKGSLAFNLQGGPPSPSLAPDGHTPFPPSPALPT THGSPQLHPLFPDPSTTMPNSTAPHPVTMYPHRNLSQETHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only