

Recombinant Human SLAM Family Member 5,SLAMF5,CD84(C-6His)

Catalog No: #AP76427

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP76427-1 10ug #AP76427-2 50ug #AP76427-3 500ug #AP76427-4 1mg

Description

Product Name	Recombinant Human SLAM Family Member 5,SLAMF5,CD84(C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9UIB8
Uniprot	Q9UIB8
GeneID	8832;
Calculated MW	23.1kDa
Target Sequence	KDSEIFTVNGILGESVTFPVNIQEPRQVKIIAWTSKTSVAYVTPGDSETAPVVTVTHRNYERIHALGPNYNLVIS DLRMEDAGDYKADINTQADPYTTTKRYNLQIYRRLGKPKITQSLMASVNSTCNVTLTCSVEKEEKNVTYNWSP LGEEGNVLQIFQTPEDQELTYTCTAQNPVSNNSDSISARQLCADIAMGFRHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only