

Recombinant Rat Interleukin-12,IL-12 (C-6His)

Catalog No: #AP76432



Package Size: #AP76432-1 10ug #AP76432-2 50ug #AP76432-3 500ug #AP76432-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Rat Interleukin-12,IL-12 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Rat
Accession No.	Q9R278&Q9R103
Uniprot	Q9R278
Calculated MW	36.7&22.6kDa
Target Sequence	MWELEKDVYVVEVDWRPDAPGETVTLTCDSPPEEDITWTSQRRGVIGSGKTLTITVREFLDAGQYTCHRGG ETLSHSHLLLHKKENGIWSTEILKNFKNKTKFLKCEAPNYSGRFTCSWLVRNTDLKFNKSSSSSPESRAVTCG RASLSAEKVTLNQRDYEKYSVACQEDVTCPTAEETLPIELVVEAQQQNKYENYSTSFFIRDIIPDPKLNQVK PLKNSQVEVSWEYPDSWSTPHSYFSLKFFVRIQRKKEKTKETEEECNQKGAFLVEKTSAEVQCKGANICVQA QDRYNNSSCSKWTCVPCRGRSHHHHHH&RVIPVSGPAKCLNQSQNLLKTTDDMVRTARELKHYSCTAGDI DHEDITRDKTSTLEACLPLELHKNESCLATKETSSIIIRGSLPPQKTSMMTLCLGSIYEDLKMYQSEFQAINAA LQSHNHQQITLDRNMLMAIDELMRSLNHSGETLHQKAPMGEADPYRVKMKLCILLHAFSTRVMTINRVMNYLS SSHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only