

Recombinant Human B-cell Maturation Protein,BCMA,TNFRSF17 (C-Fc)



Catalog No: #AP76441

Orders: order@signalwayantibody.com

Package Size: #AP76441-1 10ug #AP76441-2 50ug #AP76441-3 500ug #AP76441-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human B-cell Maturation Protein,BCMA,TNFRSF17 (C-Fc)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q02223
Uniprot	Q02223
GeneID	608;
Calculated MW	32.8kDa
Target Sequence	MLQMAGQCSQNEYFDSLLHACIPCQLRCSNTPLTCQRYCNASVTNSVKGTNADIEGRMDEPKSCDKTHTC PPCPAPELLGGPSVFLFPPKPKDTLMISRTPEVTCVVVDVSHEDPEVKFNWYVDGVEVHNAKTKPREEQYNS TYRVVSVLTVLHQDWLNGKEYKCKVSNKALPAPIEKTISKAKGQPREPQVYTLPPSREEMTKNQVSLTCLVKG FYPSDIAVEWESNGQPENNYKTTTPVLDSGDSFFLYSKLTVDKSRWQQGNVFCSSVMHEALHNHYTQKSLSL SPGK
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS,pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only