

Recombinant Human Renin (C-10His)

Catalog No: #AP76442



Package Size: #AP76442-1 10ug #AP76442-2 50ug #AP76442-3 500ug #AP76442-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Renin (C-10His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P00797
Uniprot	P00797
GeneID	5972;
Calculated MW	44kDa
Target Sequence	LPTDTTTFKRIFLKRMPHSIRESLKERGVDMARLGPEWSQPMKRLTLGNTTSSVILTNYMDTQYYGEIGITPPQ TFKVVFDTGSSNVVWPSSKCSRLYTACVYHKLFDASDSSSYKHNGTELTLYSTGTVSGFLSQDIITVGGITVT QMFGEVTEMPALPFMLAEFDG VVGMGFIEQAIGRVTPIFDNIISQGV LKEDVFSFYNRDSENSQSLGGQIVLG GSDPQHYEGNFHYINLIKTGVWQIQMKGVSVGSSTLLCEDGCLALVDTGASYISGSTSSIEKLMEALGAKKRLF DYVVKCNEGPTLPDISFHLGGKEYTLTSADYVFQESYSSKKLCTLAIHAMDIPPPTGPTWALGATFIRKFYTEFD RRNNRIGFALARGGGGSHHHHHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only