

Recombinant Human CD226 Antigen, DNAM-1, CD226 (C-6His)

Catalog No: #AP76444

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP76444-1 10ug #AP76444-2 50ug #AP76444-3 500ug #AP76444-4 1mg

Description

Product Name	Recombinant Human CD226 Antigen, DNAM-1, CD226 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q15762
Uniprot	Q15762
GeneID	10666;
Calculated MW	26.8kDa
Target Sequence	EEVLWHTSVPF AENMSLECVYPSMGILTQVEWFKIGTQQDSIAIFSP THGMVIRKPYAERVYFLNSTMASNNM TLFFRNASEDDVGYSCSLYTPQGTWQKVIQVQSDSFEAAVPSNSHIVSEPGKNVLTTCQPQMTWPVQAV RWEKIQPRQIDLLTYCNLVHGRNFTSKFPRQIVSNCSHGRWSVIVIPDVTVSDSGLYRCYLQASAGENETFVM RLTVAEGKTDNHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only