

Recombinant Human Signal-regulatory Protein Beta-1,SIRPB1,CD172b (C-6His)

Catalog No: #AP76447

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP76447-1 10ug #AP76447-2 50ug #AP76447-3 500ug #AP76447-4 1mg

Description

Product Name	Recombinant Human Signal-regulatory Protein Beta-1,SIRPB1,CD172b (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	O00241
Uniprot	O00241
GeneID	10326;
Calculated MW	37.4kDa
Target Sequence	EDELQVIQPEKSVSVAAGESATLRCAMTSLIPVGPIMWFRGAGAGRELIYNQKEGHFPRVTTVSELTKRNLD FSISINITPADAGTYCYVKFRKGGSPDDVEFKSGAGTELSVRAKPSAPVVS GP A VRATPEHTVSFTCESHGFS RDITLKWFKNGNELSDFQTNVDPAGDSVSYSIHSTARVVL TRGDVHSQVICEIAHITLQGDPLRGTANLSEAIRV PPTLEVTQQPMRAENQANVTCQVSNFYPRGLQLTWLENGNVSRTETASTLIENKDGTYNWMSWLLVNTCAH RDDVVLTCQVEHDGQQAVSKSYALEISAHQKEHGS DITHEGSHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only