

# Recombinant Mouse Chitinase-Like Protein 3 ,YM1,ECF-L,CHI3L3 (C-6His)

Catalog No: #AP76471

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP76471-1 10ug #AP76471-2 50ug #AP76471-3 500ug #AP76471-4 1mg

## Description

Product Name	Recombinant Mouse Chitinase-Like Protein 3 ,YM1,ECF-L,CHI3L3 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	O35744
Uniprot	O35744
GeneID	12655;
Calculated MW	43.4kDa
Target Sequence	YQLMCYYTSWAKDRPIEGSFKPGNIDPCLCTHLYAFAGMQNNEITYTHEQDLRDYEALNGLKDKNTELKTLA IGGWKFGPASFSAMVSTPQNRQIFIQSVIRFLRQYNFDGLNLDWQYPGSRGSPPKDKHLFSVLVKEMRKAFE EESVEKDIPRLLLTSTGAGIIDVIKSGYKIPELSQSLDIQVMTYDLHDPKDGYTGENSPLYKSPYDIGKSADLNV DSIISYWKDHGAASEKLIVGFPAYGHTFILSDPSKTGIGAPTISTGPPGKYTESGLLAYVEVCTFLNEGATEVW DAPQEVPIYAYQGNEWWGYDNVRSFKLKAQWLKDNNLGGAVVWPLDMDDFSGSFCHQRHFPLTSTLKGDLNI HSASCKGPYVDHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only