

HNF-3 α / β / γ (Acetyl-Lys264/253/211) Polyclonal Antibody

Catalog No: #HW115

Orders: order@signalwayantibody.com

Package Size: #HW115-1 50ul #HW115-2 100ul

Support: tech@signalwayantibody.com

Description

Product Name	HNF-3 α / β / γ (Acetyl-Lys264/253/211) Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	Acetyl-HNF-3 α / β / γ (K264/253/211) Polyclonal Antibody detects endogenous levels of HNF-3 α / β / γ protein only when acetylated at K264/253/211.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the human HNF-3 α / β / γ around the acetylation site of K264/253/211.
Target Name	HNF-3 α / β / γ
Modification	Acetyl
Other Names	FOXA1; HNF3A; TCF3A; Hepatocyte nuclear factor 3-alpha; HNF-3-alpha; HNF-3A; Forkhead box protein A1; Transcription factor 3A; TCF-3A; FOXA2; HNF3B; TCF3B; Hepatocyte nuclear factor 3-beta; HNF-3-beta; HNF-3B;
Accession No.	Swiss-Prot: P55317NCBI Gene ID: 3169
Uniprot	P55317
GeneID	3169;
Target Species	human
SDS-PAGE MW	50kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western blotting: 1/500 - 1/2000


ELISA: 1/20000

Not yet tested in other applications

Images



Western Blot analysis of JK cells using Acetyl-HNF-3α/β/γ (K264/253/211) Polyclonal Antibody



Western Blot analysis of Jurkat cells using Acetyl-HNF-3α/β/γ (K264/253/211) Polyclonal Antibody

Note: This product is for in vitro research use only