

Recombinant Escherichia coli Protein FimH(fimH)

Catalog No: #AP76639



Package Size: #AP76639-1 20ug #AP76639-2 100ug #AP76639-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Escherichia coli Protein FimH(fimH)
Brief Description	Recombinant Protein
Host Species	Yeast
Purification	Greater than 90% as determined by SDS-PAGE.
Immunogen Description	Expression Region:22-300aaSequence Info:Full Length
Accession No.	P08191
Uniprot	P08191
GeneID	948847;
Calculated MW	29.1 kDa
Tag Info	NO-tagged
Target Sequence	FACKTANGTAIPIGGGSANVYVNLAPVVNVGQNLVVDLSTQIFCHNDYPETITDYVTLQRGSAYGGVLSNFSGT VKYSGSSYPFPTTSETPRVVYNSRTDKPWPVALYLTPVSSAGGVAIKAGSLIAVLIRQTNNYNSDDFQFVWNI YANNDVVVPTGGCDVSARDVTVTLDPYPGSVPIPLTVYCAKSQNLGYLSGTTADAGNSIFTNTASFSPAQGV GVQLTRNGTIIPANNTVSLGAVGTSVAVSLGLTANYARTGGQVTAGNVQSIIGVTFVYQ
Formulation	Tris-based buffer50% glycerol
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C,-80°C. The shelf life of lyophilized form is 12 months at -20°C,-80°C.Notes:Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Involved in regulation of length and mediation of adhesion of type 1 fimbriae (but not necessary for the production of fimbriae). Adhesin responsible for the binding to D-mannose. It is laterally positioned at intervals in the structure of the type 1 fimbriae. In order to integrate FimH in the fimbriae FimF and FimG are needed.

References

"Highly accurate genome sequences of Escherichia coli K-12 strains MG1655 and W3110."Hayashi K., Morooka N., Yamamoto Y., Fujita K., Isono K., Choi S., Ohtsubo E., Baba T., Wanner B.L., Mori H., Horiuchi T.Mol. Syst. Biol. 2:E1-E5(2006).Research Topic:Microbiology

Note: This product is for in vitro research use only