

HSP 90 (Acetyl-Lys292/284) Polyclonal Antibody

Catalog No: #HW118



Package Size: #HW118-1 50ul #HW118-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	HSP 90 (Acetyl-Lys292/284) Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	Acetyl-HSP 90 (K292/284) Polyclonal Antibody detects endogenous levels of HSP 90 protein only when acetylated at K292/284.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the human HSP 90 around the acetylation site of K292/284.
Target Name	HSP 90
Modification	Acetyl
Other Names	HSP90AA1; HSP90A; HSPC1; HSPCA; Heat shock protein HSP 90-alpha; Heat shock 86 kDa; HSP 86; HSP86; Renal carcinoma antigen NY-REN-38; HSP90AB1; HSP90B; HSPC2; HSPCB; Heat shock protein HSP 90-beta; HSP 90; Heat shock 84 kDa; HSP 84; HSP84
Accession No.	Swiss-Prot: P07900NCBI Gene ID: 3326
Uniprot	P07900
GeneID	3320;
Target Species	human
SDS-PAGE MW	85kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western blotting: 1/500 - 1/2000


ELISA: 1/20000

Not yet tested in other applications

Images



Western Blot analysis of K562 cells using Acetyl-HSP 90 (K292/284) Polyclonal Antibody



Western Blot analysis of K562 cells using Acetyl-HSP 90 (K292/284) Polyclonal Antibody

Note: This product is for in vitro research use only